


DEPT135-mmcd-clem D2O /opt/topspin216 mfjofre {3 C7 - 331}

179.17

53.88
50.56
48.08
47.27
46.42
45.80
45.74
45.60
44.77
43.92
43.12
40.63
37.31



NAME cq_00115-100mm-d2o-10-29-10
EXPNO 5
PROCNO 1
Date_ 20101029
Time 12.40
INSTRUM spect
PROBHD 5 mm CPTXO 13C
PULPROG dept135
TD 65536
SOLVENT D2O
NS 256
DS 4
SWH 29761.904 Hz
FIDRES 0.454131 Hz
AQ 1.1010548 sec
RG 203
DW 16.800 usec
DE 30.00 usec
TE 298.0 K
CNST2 145.0000000
D1 2.00000000 sec
D2 0.00344828 sec
D12 0.00002000 sec
TD0 1

=====
CHANNEL f1
NUC1 13C
P1 10.38 usec
P2 20.76 usec
PL1 3.00 dB
PL1W 36.25705338 W
SFO1 125.6974365 MHz

=====
CHANNEL f2
CPDPRG2 waltz16
NUC2 1H
P3 15.70 usec
P4 31.40 usec
PCPD2 100.00 usec
PL2 1.25 dB
PL12 17.33 dB
PL2W 19.96668053 W
PL12W 0.49238619 W
SFO2 499.8419994 MHz
SI 32768
SF 125.6848680 MHz
WDW EM
SSB 0
LB 1.00 Hz
GB 0
PC 1.40

