



HSQC-mmcd-Ravi D2O /opt/topspin216 mfjofre {3 D7 - 343}



NAME cq_00066-100mm-d2o-10-29-10
EXPNO 6
PROCNO 1
Date_ 20101031
Time 5.18
INSTRUM spect
PROBHD 5 mm CPTXO 13C
PULPROG hsqcetgp
TD 1024
SOLVENT D2O
NS 8
DS 16
SWH 6493.506 Hz
FIDRES 6.341315 Hz
AQ 0.0788980 sec
RG 203
DW 77.000 usec
DE 30.00 usec
TE 298.0 K
CNST2 145.0000000
D0 0.00000300 sec
D1 2.00000000 sec
D4 0.00172414 sec
D11 0.03000000 sec
D13 0.00000400 sec
D16 0.00020000 sec
INO 0.00002650 sec
ZGPTNS
===== CHANNEL f1 =====
NUC1 1H
P1 15.00 usec
P2 30.00 usec
P28 0.10 usec
PL1 0.75 dB
PL1W 22.40298462 W
SFO1 499.8423492 MHz
===== CHANNEL f2 =====
CPDPRG2 garp
NUC2 13C
P3 10.00 usec
P4 20.00 usec
PCPD2 70.00 usec
PL2 3.00 dB
PL12 19.90 dB
PL2W 36.25705338 W
PL12W 0.74027407 W
SFO2 125.6936659 MHz
===== GRADIENT CHANNEL =====
GPNAM1 SINE.100
GPNAM2 SINE.100
GPZ1 80.00 %
GPZ2 20.10 %
P16 1000.00 usec
ND0 2
TD 256
SFO1 125.6937 MHz
FIDRES 73.648636 Hz
SW 150.000 ppm
FnMODE Echo-Antiecho
SI 1024
SF 499.8400000 MHz
WDW QSINE
SSB 2
LB 0.00 Hz
GB 0
PC 1.40
SI 1024
MC2 echo-antiecho
SF 125.6848680 MHz
WDW QSINE
SSB 2
LB 0.00 Hz
GB 0