


DEPT135-mmcd-clem D2O /opt/topspin216 mfjofre {3 D7 - 343}

—178.28

38.38  
36.36  
33.46  
33.41  
33.19  
28.45



NAME cq\_00066-100mm-d2o-10-29-10  
EXPNO 5  
PROCNO 1  
Date\_ 20101031  
Time 5.17  
INSTRUM spect  
PROBHD 5 mm CPTXO 13C  
PULPROG dept135  
TD 65536  
SOLVENT D2O  
NS 256  
DS 4  
SWH 29761.904 Hz  
FIDRES 0.454131 Hz  
AQ 1.1010548 sec  
RG 203  
DW 16.800 usec  
DE 30.00 usec  
TE 298.0 K  
CNST2 145.0000000  
D1 2.00000000 sec  
D2 0.00344828 sec  
D12 0.00002000 sec  
TD0 1

===== CHANNEL f1 =====  
NUC1 13C  
P1 10.38 usec  
P2 20.76 usec  
PL1 3.00 dB  
PL1W 36.25705338 W  
SFO1 125.6974365 MHz

===== CHANNEL f2 =====  
CPDPRG2 waltz16  
NUC2 1H  
P3 15.70 usec  
P4 31.40 usec  
PCPD2 100.00 usec  
PL2 1.25 dB  
PL12 17.33 dB  
PL2W 19.96668053 W  
PL12W 0.49238619 W  
SFO2 499.8419994 MHz  
SI 32768  
SF 125.6848680 MHz  
WDW EM  
SSB 0  
LB 1.00 Hz  
GB 0  
PC 1.40

